



Propranolol [HRP] (DAGA-358H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Propranolol is a non-selective beta-adrenergic receptor blocking agent. It has no other autonomic nervous system activity. Propranolol is a competitive antagonist which specifically competes with beta-adrenergic receptor stimulating agents for available beta-receptor sites. When access to beta-adrenergic receptor sites is blocked by propranolol, the chronotropic, inotropic, and vasodilator responses to beta-adrenergic stimulation are decreased proportionately.
Conjugate	HRP
Applications	ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	1 Isopropylamino 3 1 Naphthyloxy3Propanol Hydrochloride; Inderal; Inderol; Propranolol HCl; Propranolol Hydrochloride
-----------------	---