



Propranolol [KLH] (DAGA-358K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|--|
| Antigen Description | Propranolol is a non-selective beta-adrenergic receptor blocking agent. It has no other autonomic nervous system activity. Propranolol is a competitive antagonist which specifically competes with beta-adrenergic receptor stimulating agents for available beta-receptor sites. When access to beta-adrenergic receptor sites is blocked by propranolol, the chronotropic, inotropic, and vasodilator responses to beta-adrenergic stimulation are decreased proportionately. |
| Conjugate | KLH |
| Applications | Immunogen |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 1 mg |
| Buffer | 0.01M pH7.4 PBS |
| Preservative | None |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

BACKGROUND

| | |
|-----------------|---|
| Keywords | 1 Isopropylamino 3 1 Naphthyloxy3Propanol Hydrochloride; Inderal; Inderol; Propranolol HCl; Propranolol Hydrochloride |
|-----------------|---|