



Recombinant Protein Ailuropoda melanoleuca (Giant panda) RNF182, E.coli (DAGA-3606)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview

The following tags are available.

N-terminal His-tagged and C-terminal Myc-tagged
 N-terminal His-SUMO-tagged and C-terminal Myc-tagged
 N-terminal His-B2M-tagged and C-terminal Myc-tagged
 N-terminal His-GST-tagged and C-terminal Myc-tagged
 N-terminal His-tagged
 N-terminal His-SUMO-tagged
 N-terminal His-B2M-tagged
 N-terminal His-GST-tagged
 N-terminal GST-tagged
 Tag-Free

Antigen Description

Recombinant 6xHis tagged Ailuropoda melanoleuca (Giant panda) RNF182 (UniProtKB No. D2H788) was expressed.

Species

Ailuropoda melanoleuca

Purity

>85% (SDS-PAGE)

Conjugate

His

Format

Liquid

Concentration

Batch dependent - please inquire should you have specific requirements.

Size

1 mg

Buffer

Tris-based buffer, 50% glycerol

Storage

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to

one week.
