



## Theophylline [KLH] (DAGA-360K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Theophylline is a colorless crystalline alkaloid, derived from tea leaves or made synthetically, used in medicine especially as a bronchial dilator.
<b>Conjugate</b>	KLH
<b>Applications</b>	Immunogen
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

<b>Keywords</b>	Dimethylxanthine; theophylline; accurbron; Acet-theocin; aerolate; Aerolate III; armophylline; austyn; bronkodyl; bronkodylsr; choledylsa; Pseudotheophylline
-----------------	---