



Recombinant Protein Ailurus fulgens (Lesser panda) (Red panda) MT-ND2, E.coli (DAGA-3610)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview

The following tags are available.

N-terminal His-tagged and C-terminal Myc-tagged
 N-terminal His-SUMO-tagged and C-terminal Myc-tagged
 N-terminal His-B2M-tagged and C-terminal Myc-tagged
 N-terminal His-GST-tagged and C-terminal Myc-tagged
 N-terminal His-tagged
 N-terminal His-SUMO-tagged
 N-terminal His-B2M-tagged
 N-terminal His-GST-tagged
 N-terminal GST-tagged
 Tag-Free

Antigen Description	Recombinant 6xHis tagged Ailurus fulgens (Lesser panda) (Red panda) MT-ND2 (UniProtKB No. Q534E0) was expressed.
----------------------------	--

Species	Ailurus fulgens
----------------	-----------------

Purity	>85% (SDS-PAGE)
---------------	-----------------

Conjugate	His
------------------	-----

Format	Liquid
---------------	--------

Concentration	Batch dependent - please inquire should you have specific requirements.
----------------------	---

Size	1 mg
-------------	------

Buffer	Tris-based buffer, 50% glycerol
---------------	---------------------------------

Storage

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
