



Valproic acid [HRP] (DAGA-361H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Valproic acid (VPA, valproate), an acidic chemical compound, has found clinical use as an anticonvulsant and mood-stabilizing drug, primarily in the treatment of epilepsy, bipolar disorder and prevention of migraine headaches. VPA is a liquid at room temp
Conjugate	HRP
Applications	ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	VPA; Valproic acid
-----------------	--------------------