



## Valproic acid [KLH] (DAGA-361K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Antigen Description	Valproic acid (VPA, valproate), an acidic chemical compound, has found clinical use as an anticonvulsant and mood-stabilizing drug, primarily in the treatment of epilepsy, bipolar disorder and prevention of migraine headaches.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

Keywords	VPA; Valproic acid
----------	--------------------