



# Recombinant Protein *Ajellomyces capsulatus* (strain NAm1 / WU24) VMA21, Baculovirus (DAGA-3627)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |                                                                                                                                              |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Overview</b>    | <p>The following tags are available.</p> <p>N-terminal His-tagged and C-terminal Myc-tagged</p> <p>N-terminal His-tagged</p> <p>Tag-Free</p> |
| <b>Antigen Description</b> | Recombinant 6xHis tagged <i>Ajellomyces capsulatus</i> (strain NAm1 / WU24) VMA21 (UniProtKB No. A6R3V7) was expressed.                      |
| <b>Species</b>             | <i>Ajellomyces capsulatus</i>                                                                                                                |
| <b>Purity</b>              | >85% (SDS-PAGE)                                                                                                                              |
| <b>Conjugate</b>           | His                                                                                                                                          |
| <b>Format</b>              | Liquid                                                                                                                                       |
| <b>Concentration</b>       | Batch dependent - please inquire should you have specific requirements.                                                                      |
| <b>Size</b>                | 1 mg                                                                                                                                         |
| <b>Buffer</b>              | Tris-based buffer, 50% glycerol                                                                                                              |
| <b>Storage</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                          |