



# Recombinant Protein *Ajellomyces capsulatus* (strain NAm1 / WU24) CBR1, Mammalian cell (DAGA-3628)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	<p>The following tags are available.</p> <p>N-terminal His-tagged and C-terminal Myc-tagged</p> <p>N-terminal His-tagged</p> <p>Tag-Free</p>
<b>Antigen Description</b>	Recombinant 6xHis tagged <i>Ajellomyces capsulatus</i> (strain NAm1 / WU24) CBR1 (UniProtKB No. A6R2K7) was expressed.
<b>Species</b>	<i>Ajellomyces capsulatus</i>
<b>Purity</b>	>85% (SDS-PAGE)
<b>Conjugate</b>	His
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	Tris-based buffer, 50% glycerol
<b>Storage</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.