



## Warfarin [KLH] (DAGA-362K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Antigen Description

Warfarin (also known under the brand names Coumadin, Jantoven, Marevan, Lawarin, Waran, and Warfant) is an anticoagulant normally used in the prevention of thrombosis and thromboembolism, the formation of blood clots in the blood vessels and their migration elsewhere in the body respectively. It was initially introduced in 1948 as a pesticide against rats and mice and is still used for this purpose, although more potent poisons such as brodifacoum have since been developed. In the early 1950s warfarin was found to be effective and relatively safe for preventing thrombosis and embolism (abnormal formation and migration of blood clots) in many disorders. It was approved for use as a medication in 1954 and has remained popular ever since; warfarin is the most widely prescribed oral anticoagulant drug in North America.

Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

#### Keywords

Coumadin; Jantoven; Marevan; Lawarin; Waran; and Warfant; Warfarin