



Recombinant Protein *Akkermansia muciniphila* (strain ATCC BAA-835 / Muc) Amuc_0430, Baculovirus (DAGA-3636)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The following tags are available.
	N-terminal His-tagged and C-terminal Myc-tagged
	N-terminal His-tagged
	Tag-Free
	The tag type will be determined during production process. If you have specified tag type, please tell us and we will develop the specified tag preferentially. If we fail in developing the specified tag, we will send other tag type directly without negotiation.
Antigen Description	Recombinant 6xHis tagged <i>Akkermansia muciniphila</i> (strain ATCC BAA-835 / Muc) Amuc_0430 (UniProtKB No. B2UNF6) was expressed.
Species	<i>Akkermansia muciniphila</i>
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer, 50% glycerol
Storage	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.