



Recombinant Protein *Alcaligenes xylosoxydans xylosoxydans dh1C*, Mammalian cell (DAGA-3654)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The following tags are available. N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free The tag type will be determined during production process. If you have specified tag type, please tell us and we will develop the specified tag preferentially. If we fail in developing the specified tag, we will send other tag type directly without negotiation.
Antigen Description	Recombinant 6xHis tagged <i>Alcaligenes xylosoxydans xylosoxydans dh1C</i> (UniProtKB No. Q59168) was expressed.
Species	<i>Alcaligenes xylosoxydans xylosoxydans</i>
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer, 50% glycerol
Storage	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.