



Recombinant Protein *Alicyclobacillus acidocaldarius* subsp. *acidocaldarius* (strain ATCC 27009) kdpA, Baculovirus (DAGA-3673)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	<p>The following tags are available.</p> <p>N-terminal His-tagged and C-terminal Myc-tagged</p> <p>N-terminal His-tagged</p> <p>Tag-Free</p>
Antigen Description	Recombinant 6xHis tagged <i>Alicyclobacillus acidocaldarius</i> subsp. <i>acidocaldarius</i> (strain ATCC 27009) kdpA (UniProtKB No. Q9XE11) was expressed.
Species	<i>Alicyclobacillus acidocaldarius</i> subsp. <i>acidocaldarius</i>
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer, 50% glycerol
Storage	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.