



Recombinant Protein Allochrocatium vinosum (strain ATCC 1789) petC, Yeast (DAGA-3711)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The following tags are available. N-terminal His-tagged Tag-Free
Antigen Description	Recombinant 6xHis tagged Allochrocatium vinosum (strain ATCC 1789) petC (UniProtKB No. O31216) was expressed.
Species	Allochrocatium vinosum
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer,50% glycerol
Storage	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.