



Recombinant Protein Allpahuayo mammarenavirus L, Mammalian cell (DAGA- 3723)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview

The following tags are available.

N-terminal His-tagged and C-terminal Myc-tagged

N-terminal His-tagged

Tag-Free

The tag type will be determined during production process. If you have specified tag type, please tell us and we will develop the specified tag preferentially. If we fail in developing the specified tag, we will send other tag type directly without negotiation.

Antigen Description

Recombinant 6xHis tagged Allpahuayo mammarenavirus L (UniProtKB No. Q6XQI9) was expressed.

Species

Allpahuayo mammarenavirus

Purity

>85% (SDS-PAGE)

Conjugate

His

Format

Liquid

Concentration

Batch dependent - please inquire should you have specific requirements.

Size

1 mg

Buffer

Tris-based buffer, 50% glycerol

Storage

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.