



Recombinant Protein *Alteromonas mediterranea* uppP, E.coli (DAGA-3753)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview

The following tags are available.

N-terminal His-tagged and C-terminal Myc-tagged
 N-terminal His-SUMO-tagged and C-terminal Myc-tagged
 N-terminal His-B2M-tagged and C-terminal Myc-tagged
 N-terminal His-GST-tagged and C-terminal Myc-tagged
 N-terminal His-tagged
 N-terminal His-SUMO-tagged
 N-terminal His-B2M-tagged
 N-terminal His-GST-tagged
 N-terminal GST-tagged
 Tag-Free

Antigen Description	Recombinant 6xHis tagged <i>Alteromonas mediterranea</i> uppP (UniProtKB No. B4RY28) was expressed.
----------------------------	---

Species	<i>Alteromonas mediterranea</i>
----------------	---------------------------------

Purity	>85% (SDS-PAGE)
---------------	-----------------

Conjugate	His
------------------	-----

Format	Liquid
---------------	--------

Concentration	Batch dependent - please inquire should you have specific requirements.
----------------------	---

Size	1 mg
-------------	------

Buffer	Tris-based buffer, 50% glycerol
---------------	---------------------------------

Storage	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to
----------------	---

one week.
