



Recombinant Protein *Alteromonas naphthalenivorans* ftsW, E.coli (DAGA-3757)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview

The following tags are available.

- N-terminal His-tagged and C-terminal Myc-tagged
- N-terminal His-SUMO-tagged and C-terminal Myc-tagged
- N-terminal His-B2M-tagged and C-terminal Myc-tagged
- N-terminal His-GST-tagged and C-terminal Myc-tagged
- N-terminal His-tagged
- N-terminal His-SUMO-tagged
- N-terminal His-B2M-tagged
- N-terminal His-GST-tagged
- N-terminal GST-tagged
- Tag-Free

Antigen Description	Recombinant 6xHis tagged <i>Alteromonas naphthalenivorans</i> ftsW (UniProtKB No. F5ZBR5) was expressed.
Species	<i>Alteromonas naphthalenivorans</i>
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer, 50% glycerol
Storage	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to

one week.
