



Recombinant Protein Amphidinium operculatum psaA, E.coli (DAGA-3809)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview

The following tags are available.

N-terminal His-tagged and C-terminal Myc-tagged
 N-terminal His-SUMO-tagged and C-terminal Myc-tagged
 N-terminal His-B2M-tagged and C-terminal Myc-tagged
 N-terminal His-GST-tagged and C-terminal Myc-tagged
 N-terminal His-tagged
 N-terminal His-SUMO-tagged
 N-terminal His-B2M-tagged
 N-terminal His-GST-tagged
 N-terminal GST-tagged
 Tag-Free

Antigen Description Recombinant 6xHis tagged Amphidinium operculatum psaA (UniProtKB No. Q9MTQ4) was expressed.

Species Amphidinium operculatum

Purity >85% (SDS-PAGE)

Conjugate His

Format Liquid

Concentration Batch dependent - please inquire should you have specific requirements.

Size 1 mg

Buffer Tris-based buffer, 50% glycerol

Storage Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to

one week.
