



Recombinant Protein *Anabaena variabilis* psaJ, *E.coli* (DAGA-3834)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview

The following tags are available.

N-terminal His-tagged and C-terminal Myc-tagged
 N-terminal His-SUMO-tagged and C-terminal Myc-tagged
 N-terminal His-B2M-tagged and C-terminal Myc-tagged
 N-terminal His-GST-tagged and C-terminal Myc-tagged
 N-terminal His-tagged
 N-terminal His-SUMO-tagged
 N-terminal His-B2M-tagged
 N-terminal His-GST-tagged
 N-terminal GST-tagged
 Tag-Free

Antigen Description	Recombinant 6xHis tagged <i>Anabaena variabilis</i> psaJ (UniProtKB No. P23080) was expressed.
Species	<i>Anabaena variabilis</i>
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer, 50% glycerol
Storage	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.