



Recombinant Protein Anaeromyxobacter dehalogenans (strain 2CP-1 / ATCC BAA-258) atpB, E.coli (DAGA-3844)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	<p>The following tags are available.</p> <p>N-terminal His-tagged and C-terminal Myc-tagged</p> <p>N-terminal His-SUMO-tagged and C-terminal Myc-tagged</p> <p>N-terminal His-B2M-tagged and C-terminal Myc-tagged</p> <p>N-terminal His-GST-tagged and C-terminal Myc-tagged</p> <p>N-terminal His-tagged</p> <p>N-terminal His-SUMO-tagged</p> <p>N-terminal His-B2M-tagged</p> <p>N-terminal His-GST-tagged</p> <p>N-terminal GST-tagged</p> <p>Tag-Free</p>
Antigen Description	Recombinant 6xHis tagged Anaeromyxobacter dehalogenans (strain 2CP-1 / ATCC BAA-258) atpB (UniProtKB No. B8JCU2) was expressed.
Species	Anaeromyxobacter dehalogenans
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer, 50% glycerol

Storage

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
