



Recombinant Protein Anaeromyxobacter dehalogenans (strain 2CP-1 / ATCC BAA-258) atpB, Mammalian cell (DAGA-3845)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The following tags are available. N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free
Antigen Description	Recombinant 6xHis tagged Anaeromyxobacter dehalogenans (strain 2CP-1 / ATCC BAA-258) atpB (UniProtKB No. B8JCU2) was expressed.
Species	Anaeromyxobacter dehalogenans
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer, 50% glycerol
Storage	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.