



# Recombinant Protein Anaeromyxobacter dehalogenans (strain 2CP-C) nhaA, E.coli (DAGA-3847)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

### Product Overview

The following tags are available.

N-terminal His-tagged and C-terminal Myc-tagged  
 N-terminal His-SUMO-tagged and C-terminal Myc-tagged  
 N-terminal His-B2M-tagged and C-terminal Myc-tagged  
 N-terminal His-GST-tagged and C-terminal Myc-tagged  
 N-terminal His-tagged  
 N-terminal His-SUMO-tagged  
 N-terminal His-B2M-tagged  
 N-terminal His-GST-tagged  
 N-terminal GST-tagged  
 Tag-Free

<b>Antigen Description</b>	Recombinant 6xHis tagged Anaeromyxobacter dehalogenans (strain 2CP-C) nhaA (UniProtKB No. Q2IE76) was expressed.
<b>Species</b>	Anaeromyxobacter dehalogenans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Conjugate</b>	His
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	Tris-based buffer, 50% glycerol

**Storage**

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

---