



Recombinant Protein Anaeromyxobacter sp. (strain K) rplV, Baculovirus (DAGA-3855)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|---|
| Product Overview | The following tags are available. |
| | N-terminal His-tagged and C-terminal Myc-tagged |
| | N-terminal His-tagged |
| | Tag-Free |
| | The tag type will be determined during production process. If you have specified tag type, please tell us and we will develop the specified tag preferentially. If we fail in developing the specified tag, we will send other tag type directly without negotiation. |
| Antigen Description | Recombinant 6xHis tagged Anaeromyxobacter sp. (strain K) rplV (UniProtKB No. B4UBA4) was expressed. |
| Species | Anaeromyxobacter |
| Purity | >85% (SDS-PAGE) |
| Conjugate | His |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 1 mg |
| Buffer | Tris-based buffer, 50% glycerol |
| Storage | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |