



Recombinant Anoxybacillus flavithermus Aflv_2299, Mammalian cell (DAGA-3941)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|---|
| Product Overview | <p>The following tags are available.</p> <p>N-terminal His-tagged and C-terminal Myc-tagged</p> <p>N-terminal His-tagged</p> <p>Tag-Free</p> |
| Antigen Description | Recombinant 6xHis tagged Anoxybacillus flavithermus Aflv_2299 (UniProtKB No. B7GF17) was expressed. |
| Purity | >85% (SDS-PAGE) |
| Conjugate | His |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 1 mg |
| Buffer | Tris-based buffer, 50% glycerol |
| Storage | <p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p> |