



Recombinant Antonospora locustae Isoleucine--tRNA ligase , partial, E. Coli (DAGA-3984)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview

The following tags are available.

N-terminal His-tagged and C-terminal Myc-tagged

N-terminal His-SUMO-tagged and C-terminal Myc-tagged

N-terminal His-B2M-tagged and C-terminal Myc-tagged

N-terminal His-GST-tagged and C-terminal Myc-tagged

N-terminal His-tagged

N-terminal His-SUMO-tagged

N-terminal His-B2M-tagged

N-terminal His-GST-tagged

N-terminal GST-tagged

Tag-Free

The tag type will be determined during production process. If you have specified tag type, please tell us and we will develop the specified tag preferentially. If we fail in developing the specified tag, we will send other tag type directly without negotiation.

Antigen Description	Recombinant 6xHis tagged Antonospora locustae Isoleucine--tRNA ligase (UniProtKB No. Q27707) was expressed.
----------------------------	---

Purity	>85% (SDS-PAGE)
---------------	-----------------

Conjugate	His
------------------	-----

Format	Liquid
---------------	--------

Concentration	Batch dependent - please inquire should you have specific requirements.
----------------------	---

Size	1 mg
-------------	------

Buffer	Tris-based buffer,50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.