



Recombinant Aquareovirus A (isolate Chum salmon/Japan/CSRV/1981) S3, Baculovirus (DAGA-4054)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The following tags are available. N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free The tag type will be determined during production process. If you have specified tag type, please tell us and we will develop the specified tag preferentially. If we fail in developing the specified tag, we will send other tag type directly without negotiation.
Antigen Description	Recombinant 6xHis tagged Aquareovirus A (isolate Chum salmon/Japan/CSRV/1981) S3 (UniProtKB No. Q8VA41) was expressed.
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer, 50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.