



Recombinant *Arcobacter nitrofigilis* cowN, E.coli (DAGA-4127)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview

The following tags are available.

N-terminal His-tagged and C-terminal Myc-tagged
N-terminal His-SUMO-tagged and C-terminal Myc-tagged
N-terminal His-B2M-tagged and C-terminal Myc-tagged
N-terminal His-GST-tagged and C-terminal Myc-tagged
N-terminal His-tagged
N-terminal His-SUMO-tagged
N-terminal His-B2M-tagged
N-terminal His-GST-tagged
N-terminal GST-tagged
Tag-Free

The tag type will be determined during production process. If you have specified tag type, please tell us and we will develop the specified tag preferentially. If we fail in developing the specified tag, we will send other tag type directly without negotiation.

Antigen Description

Recombinant 6xHis tagged *Arcobacter nitrofigilis* cowN (UniProtKB No. D5V2I8) was expressed.

Purity

>85% (SDS-PAGE)

Conjugate

His

Format

Liquid

Concentration

Batch dependent - please inquire should you have specific requirements.

Size

1 mg

Buffer

Tris-based buffer, 50% glycerol

Storage

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
