



# Recombinant Arctic squirrel hepatitis virus C, Baculovirus (DAGA-4130)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

### Product Overview

The following tags are available.

N-terminal His-tagged and C-terminal Myc-tagged

N-terminal His-tagged

Tag-Free

The tag type will be determined during production process. If you have specified tag type, please tell us and we will develop the specified tag preferentially. If we fail in developing the specified tag, we will send other tag type directly without negotiation.

### Antigen Description

Recombinant 6xHis tagged Arctic squirrel hepatitis virus C (UniProtKB No. Q64896) was expressed.

### Purity

>85% (SDS-PAGE)

### Conjugate

His

### Format

Liquid

### Concentration

Batch dependent - please inquire should you have specific requirements.

### Size

1 mg

### Buffer

Tris-based buffer, 50% glycerol

### Storage

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.