



Recombinant Aromatoleum aromaticum) AZOSEA37130, Yeast (DAGA-4164)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The following tags are available. N-terminal His-tagged Tag-Free
Antigen Description	Recombinant 6xHis tagged Aromatoleum aromaticum) AZOSEA37130 (UniProtKB No. Q5NYM6) was expressed.
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer,50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.