



Recombinant Aromatoleum aromaticum) macB, E.coli (DAGA-4166)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The following tags are available.
	N-terminal His-tagged and C-terminal Myc-tagged
	N-terminal His-SUMO-tagged and C-terminal Myc-tagged
	N-terminal His-B2M-tagged and C-terminal Myc-tagged
	N-terminal His-GST-tagged and C-terminal Myc-tagged
	N-terminal His-tagged
	N-terminal His-SUMO-tagged
	N-terminal His-B2M-tagged
	N-terminal His-GST-tagged
	N-terminal GST-tagged
	Tag-Free
Antigen Description	Recombinant 6xHis tagged Aromatoleum aromaticum) macB (UniProtKB No. Q5P6D5) was expressed.
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer,50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized

form is 12 months at -20°C/-80°C.
