



Recombinant Arthroderma benhamiae CHO2, E.coli (DAGA-4200)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview

The following tags are available.

N-terminal His-tagged and C-terminal Myc-tagged
 N-terminal His-SUMO-tagged and C-terminal Myc-tagged
 N-terminal His-B2M-tagged and C-terminal Myc-tagged
 N-terminal His-GST-tagged and C-terminal Myc-tagged
 N-terminal His-tagged
 N-terminal His-SUMO-tagged
 N-terminal His-B2M-tagged
 N-terminal His-GST-tagged
 N-terminal GST-tagged
 Tag-Free

Antigen Description

Recombinant 6xHis tagged Arthroderma benhamiae CHO2 (UniProtKB No. D4AT37) was expressed.

Purity

>85% (SDS-PAGE)

Conjugate

His

Format

Liquid

Concentration

Batch dependent - please inquire should you have specific requirements.

Size

1 mg

Buffer

Tris-based buffer, 50% glycerol

Storage

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.
 Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized

form is 12 months at -20°C/-80°C.
