



Purified Adenovirus type 40 virions (DAGA-452)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Adenovirus serotypes 40 and 41 cause acute gastroenteritis primarily in children. Symptoms may include fever, diarrhea, vomiting, and abdominal pain, and last for approximately 10 days. Respiratory symptoms can occur in some individuals. The disease is usually self-limiting in immunocompetent individuals; however rare fatalities can occur in immunocompromised individuals. Asymptomatic infections are common, particularly in children.
Purity	> 95% as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	User optimized
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -80°C

BACKGROUND

Keywords Adenovirus type 40 antigen;Adenovirus protein