



Aspergillus fumigatus (DAGA-459)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Aspergillus fumigatus is a fungus of the genus Aspergillus, and is one of the most common Aspergillus species to cause disease in individuals with an immunodeficiency.
Conjugate	Unconjugated
Applications	Suitable for use in EIA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Buffer	Phosphate Buffered Saline
Preservative	None
Storage	Store at -20°C to -80°C. Avoid multiple freeze/thaw cycles

BACKGROUND

Keywords Aspergillus fumigatus;Aspergillus fumigatus antigen;Aspergillus fumigatus protein