



Recombinant Borrelia Burgdorferi p100 [His] (DAGA-472)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant Borrelia Burgdorferi p100 (p100/p83) produced in SF9 is a glycosylated, polypeptide chain having a calculated molecular mass of 77,813 Dalton. Borrelia p100 is expressed with a -6x His tag at N-terminus and purified by proprietary chromatographic techniques.
Purity	80.0% by SDS-PAGE
Conjugate	His
Molecular Weight	38.132 kD
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	5 µg
Buffer	20mM HEPES buffer pH-7.6, 250mM NaCl and 20% glycerol.
Preservative	None
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.

BACKGROUND

Keywords	Bacteria p100;Spirochaetes p100;Spirochaetales p100;Spirochaetaceae p100;Borrelia p100;Borrelia garinii p100;Borrelia burgdorferi p100;Borrelia burgdorferi p100;Borrelia burgdorferi garini p100;B. afzelii p100 protein;p100
-----------------	--