



Recombinant SARS Nucleoprotein (a.a. 1-49, 192-220) (DAGA-606)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	E.coli derived recombinant. The protein contains Nucleocapsid protein immunodominant regions
Purity	>95% by SDS-PAGE
Conjugate	TBD
Applications	Antigen in ELISA and Western blots, excellent antigen for detection of SARS with minimal specificity problems.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	50 mM Tris-HCl, 60 mM Na Cl, 50% glycerol.
Preservative	None
Storage	Upon arrival, store at -20°C

BACKGROUND

Keywords	Coronavirus Core protein;SARS;Core protein;Core protein SARS
----------	--