



# Recombinant SARS Envelope protein (DAGA-607)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Antigen Description	E.coli derived recombinant. The protein contains the E protein immunodominant regions
Purity	>95% by SDS-PAGE
Conjugate	TBD
Applications	Antigen in ELISA and Western blots, excellent antigen for detection of SARS with minimal specificity problems.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	50 mM Tris-HCl, 60 mM Na Cl, 50% glycerol.
Preservative	None
Storage	Protein is shipped at ambient temperature. Upon arrival, store at -20°C

## BACKGROUND

Keywords	Coronavirus E protein;SARS;E protein;E protein SARS
----------	---