



Recombinant SARS Spike Protein (C-terminal section) (DAGA-609)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	E.coli derived recombinant. The protein contains the C-terminal t section of the Spike protein immunodominant regions
Purity	>95% by SDS-PAGE
Conjugate	TBD
Applications	Antigen in ELISA and Western blots, excellent antigen for detection of SARS with minimal specificity problems.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	25 mM Tris-HCl, 0.4% sarcosyl, 0,25% Triton -100, 50% glycerol.
Preservative	None
Storage	Protein is shipped at ambient temperature. Upon arrival, store at -20°C

BACKGROUND

Keywords	Coronavirus Spike Mosaic Protein S;C-terminal;SARS;Spike Mosaic Protein S;C-terminal;Spike Mosaic Protein S;C-terminal;SARS
-----------------	---