



## Cat haptoglobin reference serum (DAGA-618)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Haptoglobin (Hp) is an acute-phase protein that is crucial for the elimination of free hemoglobin and the neutralization of oxidative damage.
<b>Specificity</b>	Cat
<b>Conjugate</b>	Unconjugated
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Size</b>	1 ml
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Frozen -20°C

### BACKGROUND

**Keywords** Cat;Haptoglobin;Hp;Serum