



# Chicken alpha 1-acid glycoprotein reference serum (DAGA-619)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Alpha-1-acid glycoprotein (AGP) is a highly glycosylated protein with 41–43 kDa of MW. AGP is an acute phase protein in all mammals investigated. During an acute-phase response, serum level of AGP increases several fold as well as based on the type of inflammation its glycosylation pattern is different.
<b>Specificity</b>	Chicken
<b>Conjugate</b>	Unconjugated
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Size</b>	1 ml
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Frozen -20°C

## BACKGROUND

**Keywords** Chicken; Alpha 1-Acid Glycoprotein; Serum