



Chicken SAA reference serum (DAGA-620)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	SAA protein is an acute phase apolipoprotein reactant which is produced mostly by hepatocytes and under regulation of inflammatory cytokines. SAA1 (Serum amyloid A1) protein is produced mainly in the liver and circulates in low levels in the blood.
Specificity	Chicken
Conjugate	Unconjugated
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 ml
Preservative	0.1% Sodium Azide
Storage	Frozen -20°C

BACKGROUND

Keywords Chicken;SAA;Serum