



Dog alpha 1-acid glycoprotein reference serum (DAGA-623)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Alpha-1-acid glycoprotein (AGP) is a highly glycosylated protein with 41–43 kDa of MW. AGP is an acute phase protein in all mammals investigated. During an acute-phase response, serum level of AGP increases several fold as well as based on the type of inflammation its glycosylation pattern is different.
Specificity	Canine, Dog
Conjugate	Unconjugated
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 ml
Preservative	0.1% Sodium Azide
Storage	Frozen -20°C

BACKGROUND

Keywords Dog; Alpha 1-Acid Glycoprotein; Serum