



Dog alpha 1-antitrypsin reference serum (DAGA-624)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Alpha-1 Antitrypsin is a protease inhibitor belonging to the serpin superfamily. It is known as serum trypsin inhibitor. Alpha-1 antitrypsin is also referred to as alpha-1 proteinase inhibitor (A1PI) because it inhibits a wide variety of proteases. It protects tissues from enzymes of inflammatory cells, especially neutrophil elastase.

Specificity Canine, Dog

Conjugate Unconjugated

Format Liquid

Concentration Batch dependent - please inquire should you have specific requirements

Size 1 ml

Preservative 0.1% Sodium Azide

Storage Frozen -20°C

BACKGROUND

Keywords Dog;Alpha 1-Antitrypsin;Serum