



Human alpha 1-acid glycoprotein reference serum (DAGA-637)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Alpha-1-acid glycoprotein (α 1AGp, AGP or AAG) is an acute phase (acute phase protein) plasma alpha-globulin glycoprotein and is modulated by two polymorphic genes. It is synthesized primarily in hepatocytes and has a normal plasma concentration between 0.6-1.2 mg/mL (1-3% plasma protein). Plasma levels are affected by pregnancy, burns, certain drugs, and certain diseases, particularly HIV.
Specificity	Human
Conjugate	Unconjugated
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 ml
Preservative	0.1% Sodium Azide
Storage	Frozen -20°C

BACKGROUND

Keywords	Human;Alpha 1-Acid Glycoprotein;Serum
-----------------	---------------------------------------