



# Human cystatin C reference serum (DAGA-644)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

**Antigen Description** Cystatin C or cystatin 3 (formerly gamma trace, post-gamma-globulin, or neuroendocrine basic polypeptide), a protein encoded by the CST3 gene, is mainly used as a biomarker of kidney function. It has been studied for its role in predicting new-onset or deteriorating cardiovascular disease. It also seems to play a role in brain disorders involving amyloid (a specific type of protein deposition), such as Alzheimer's disease.

**Specificity** Human

**Conjugate** Unconjugated

**Format** Liquid

**Concentration** Batch dependent - please inquire should you have specific requirements

**Size** 1 ml

**Preservative** 0.1% Sodium Azide

**Storage** Frozen -20°C

## BACKGROUND

**Keywords** Human;Cystatin C;Serum