



Mouse cystatin C reference serum (DAGA-659)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Cystatin C or cystatin 3 (formerly gamma trace, post-gamma-globulin, or neuroendocrine basic polypeptide), a protein encoded by the CST3 gene, is mainly used as a biomarker of kidney function. It has been studied for its role in predicting new-onset or deteriorating cardiovascular disease. It also seems to play a role in brain disorders involving amyloid (a specific type of protein deposition), such as Alzheimer's disease.

Specificity Mouse

Conjugate Unconjugated

Format Liquid

Concentration Batch dependent - please inquire should you have specific requirements

Size 1 ml

Preservative 0.1% Sodium Azide

Storage Frozen -20°C

BACKGROUND

Keywords Mouse;Cystatin C;Serum