



## Mouse fibrinogen reference serum (DAGA-661)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Fibrinogen (factor I) is a glycoprotein in vertebrates that helps in the formation of blood clots. It consists of a linear array of three nodules held together by a very thin thread which is estimated to have a diameter between 8 and 15 Angstrom.
<b>Specificity</b>	Mouse
<b>Conjugate</b>	Unconjugated
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Size</b>	1 ml
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Frozen -20°C

### BACKGROUND

<b>Keywords</b>	Mouse;Fibrinogen;Serum
-----------------	------------------------