



Mouse hemoglobin reference serum (DAGA-663)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Hemoglobin is the iron-containing oxygen-transport metalloprotein in the red blood cells of all vertebrates. Hemoglobin in the blood carries oxygen from the respiratory organs (lungs or gills) to the rest of the body. There it releases the oxygen to permit aerobic respiration to provide energy to power the functions of the organism in the process called metabolism.
Specificity	Mouse
Conjugate	Unconjugated
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 ml
Preservative	0.1% Sodium Azide
Storage	Frozen -20°C

BACKGROUND

Keywords Mouse;Hemoglobin;Serum