



## Pig plasminogen reference serum (DAGA-677)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                            |                                                                                                                                                                                                                                |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Antigen Description</b> | Plasmin is an important enzyme present in blood that degrades many blood plasma proteins, including fibrin clots. The degradation of fibrin is termed fibrinolysis. In humans, the plasmin protein is encoded by the PLG gene. |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|                    |              |
|--------------------|--------------|
| <b>Specificity</b> | Pig, Porcine |
|--------------------|--------------|

|                  |              |
|------------------|--------------|
| <b>Conjugate</b> | Unconjugated |
|------------------|--------------|

|               |        |
|---------------|--------|
| <b>Format</b> | Liquid |
|---------------|--------|

|                      |                                                                        |
|----------------------|------------------------------------------------------------------------|
| <b>Concentration</b> | Batch dependent - please inquire should you have specific requirements |
|----------------------|------------------------------------------------------------------------|

|             |      |
|-------------|------|
| <b>Size</b> | 1 ml |
|-------------|------|

|                     |                   |
|---------------------|-------------------|
| <b>Preservative</b> | 0.1% Sodium Azide |
|---------------------|-------------------|

|                |              |
|----------------|--------------|
| <b>Storage</b> | Frozen -20°C |
|----------------|--------------|

### BACKGROUND

|                 |                       |
|-----------------|-----------------------|
| <b>Keywords</b> | Pig;Plasminogen;Serum |
|-----------------|-----------------------|