



Rabbit alpha 2-macroglobulin reference serum (DAGA-678)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Alpha-2-Macroglobulin (α 2M) is a large plasma protein found in the blood. It is mainly produced by the liver, and also locally synthesized by macrophages, fibroblasts, and adrenocortical cells. In humans it is encoded by the A2M gene.
Specificity	Rabbit
Conjugate	Unconjugated
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 ml
Preservative	0.1% Sodium Azide
Storage	Frozen -20°C

BACKGROUND

Keywords	Rabbit;Alpha 2-Macroglobulin;Serum
-----------------	------------------------------------