



Rabbit plasminogen reference serum (DAGA-681)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Plasmin is an important enzyme present in blood that degrades many blood plasma proteins, including fibrin clots. The degradation of fibrin is termed fibrinolysis. In humans, the plasmin protein is encoded by the PLG gene.

Specificity Rabbit

Conjugate Unconjugated

Format Liquid

Concentration Batch dependent - please inquire should you have specific requirements

Size 1 ml

Preservative 0.1% Sodium Azide

Storage Frozen -20°C

BACKGROUND

Keywords Rabbit;Plasminogen;Serum